Flowgun Ultra

FLOWLIFE

Brochure

Flowgun Ultra

Percussion Massage Gun

Table of Content

4 —	Manifest
6 —	Why Flowgun Ultra
7 —	Accessories
10 —	Flowgun Ultra - A Comprehensive Guide
12 —	Getting Started
22 —	For Muscle Soreness
26 —	For General Recovery
30 —	Before, During, and After your Workout
35 —	Contact Information
36 —	Manual
	• English • Svenska • Danska • Norsk • Suomi • Deutsch • Polski • Français • Italiano • Español • Nederlands

Manifest

In a world where limits are meant to be shattered, we fuel the fire of relentless pioneers. Our brand champions those who break the mold and embrace the extraordinary, symbolizing a life lived without compromise.

Committed to relentless innovation, we design for those who demand excellence and performance, understanding that the journey into the unknown often leads to unforgettable achievements.

We want to power dreams, aspirations, and triumphs of the audacious, sparking the passions of the unstoppable.


For those who redefine boundaries and explore new horizons, we are Flowlife. Made Possible.

Why Flowgun Ultra?

- **Variable Amplitude for Targeted Treatment** - The adjustable stroke length (8–16 mm) lets you tailor each session. Use short amplitudes to gently calm sensitive areas, or go deeper to release tension and tight fascia in large muscle groups.
- **Intelligent Programs for Guided Recovery** - Four treatment modes range from steady manual control to dynamic wave patterns. Choose the rhythm that fits your goal—warming up, winding down, or stimulating circulation.
- **Rotating Arm and Triangular Frame for Ergonomic Reach** -A rotating joint and triangular frame unlock new reach and control. Switch between single- or two-handed grips to treat hard-to-reach areas with ease and precision.
- **38 Kg of Stall Force. Zero Compromise** -With up to 38 kg stall force at 8 mm and 28 kg at 16 mm, Flowgun Ultra delivers deep, consistent pressure—even in dense muscle—without losing power.
- **OLED Display with Thoughtful Charging Feedback** - A clear OLED screen displays speed, mode, amplitude, and battery status. Smart charging indicators help prevent overcharging and adapter issues.
- **Optional Heat for More Profound Relaxation** - The Heating Head adds soothing warmth up to 50°C, enhancing recovery in stiff or poorly circulated areas—before or after your session.


Accessories

Flowgun Ultra comes with everything you need for effective, full-body treatment. Each component is selected for precision, performance, and comfort, whether you're targeting deep knots or managing everyday muscle fatigue.

Flowgun Ultra

The core device features variable amplitude (8–16 mm), four treatment modes, and a rotating arm for ergonomic reach. Its triangular frame offers versatile grip for one- or two-handed use—ideal for both light and deep treatments.

Type-C Charger with EU Adapter

Fast USB-C charging for consistent uptime. Includes an EU-compatible wall adapter and built-in safeguards against overvoltage, overheating, and incorrect power source use—protecting both your device and your safety.

Carrying Case

A structured, durable case with molded compartments to protect every component and keep everything in place when on the move.


Heating Head

Combines percussion with targeted heat up to 50°C. Activates automatically during use to increase circulation and soften stiff muscles. Ideal for use on larger muscle groups.


Trigger Cone

Designed for deep, focused release in compact areas like glutes, feet, and forearms. Ideal for targeting trigger points and relieving built-up tension.


Cup Head

Rounded and evenly pressured, perfect for joint-adjacent zones and safe full-body relief.

Fork Head

A dual-point head that stimulates two fascia points simultaneously. Perfect for fascia lines that require precision and space.


Axe Head

IT bands, shoulders, and lateral muscle groups.


Silicone Hat


A soft, impact-absorbing surface that cushions each percussion. Ideal for use on sore, tender, or sensitive muscle groups — or entry-level treatment.

Flat Head

A broad, firm-surfaced head used for general recovery across thighs, calves, and larger muscle groups.

A high-contrast, black and white close-up photograph of a person's face and neck. The image focuses on the skin's texture, including pores and a small mole on the left side of the neck. The lighting creates soft shadows, emphasizing the contours of the face and neck.

Flowgun Ultra: A Comprehensive Guide

Flowgun Ultra is Flowlife's most advanced massage device to date – designed for those who demand deep relief, configurable performance, and precision control over every recovery session. Whether you're treating muscle knots, enhancing circulation, or accelerating post-workout recovery, this guide will help you master its full potential.

Getting Started

1. Getting Started

Before you begin, ensure your device is fully charged using the included USB-C cable. The charging indicator will pulse during charging and remain solid when the battery is fully charged. A complete charge takes approximately 3 hours and provides up to 3 hours of use, depending on the intensity and pressure settings.

Avoid using standard 5V chargers – only use the adapter provided.

Preparing the Device

- Attach your preferred head securely by inserting it straight into the device port.
- If using the Heating Head, it will activate automatically when the device percussion is on.

2. Understanding the Modes

Each mode on Flowgun Ultra is specifically designed for achieving distinct therapeutic outcomes. Switching modes allows you to tailor your treatment to your body's needs, whether you're activating muscles pre-workout or recovering post-session.

M1 – Manual Mode

This is the default setting for precise control. Adjust intensity in 6 levels (from 1500 to 2500 RPM) depending on your preference.

Use this mode when you want to control pressure and speed manually – ideal for specific sore spots or lighter self-massage.

M2 – Rising Mode

This mode cycles gradually upward through 6 speed levels (1500–2500 RPM), spending 5 seconds at each.

Best For

Muscle activation before training. It stimulates circulation, preps the nervous system, and primes muscles without overstimulation.

M3 – Falling Mode

The inverse of M2, this mode begins high and drops in stages. It starts at 2500 RPM and slows to 1500 RPM.

Best For

Post-workout cooldowns. It gradually winds down the muscle tone and promotes parasympathetic activation.

M4 – Rolling Wave Mode

A sinusoidal speed pattern that rapidly increases and decreases between 1700 and 2300 RPM.

Best For

Comprehensive muscle stimulation. Excellent for stubborn soreness, deep tension, or general fatigue when you're unsure where to start.

To Change Mode

Hold down the “M” button to enter mode-switching. Use short presses to cycle. The OLED screen shows your current mode and gear.

3. Adjusting Stroke Length

Every body is different — and every muscle needs a different approach. Flowgun Ultra features an adjustable stroke length mechanism (8–16 mm) that lets you control how deeply each percussion penetrates. This isn't just about intensity. It's about matching depth to muscle density, sensitivity, and your recovery goals.

What is Stroke Length?


Stroke length refers to the vertical travel distance of the massage head during each percussion. A short stroke delivers fast, surface-level stimulation. A long stroke reaches deeper into the muscle tissue, fascia, and trigger points — ideal for larger or more stubborn areas.

How Do I Adjust It?

Rotate the adjustment dial on the side of the elbow arm. This controls the amplitude of the stroke in real-time, with millimeter precision between 8 and 16 mm. You'll see clear markings along the dial to guide you.

You can set it anywhere along the spectrum — full or half increments — depending on how your body responds.

What Do the Numbers Mean?

Each number corresponds to the millimeter depth of each stroke:

8–10 mm (Short stroke)

Fast and shallow. Best for treating areas with thin or compact layers of muscle or fascia. Ideal for bony or angular regions—such as elbows, forearms, or shins—where soft tissue coverage is minimal. This setting is also suitable for warming up, boosting circulation, or working around sore or inflamed zones without overstimulating the tissue.

11–13 mm (Medium stroke)

A balanced depth for moderate tissue layers. Effective across a wide range of muscle groups, including arms, shoulders, and hamstrings. This is your go-to amplitude for general recovery—deep enough to stimulate without overwhelming, ideal when consistent pressure across medium-density tissue is desired.

14–16 mm (Long stroke)

Deep and deliberate. Best suited for thick muscle groups like the glutes, quads, and posterior chain, where more amplitude is needed to reach deeper fascial layers. These settings help release chronic tightness, improve mobility, and flush large muscle areas after intensive training or long periods of physical strain.

What happens in the body?

Shorter strokes (8–10 mm)

Activate superficial sensory receptors and stimulate circulation in surface-level tissue. This improves fluid exchange and neuromuscular readiness without stressing the underlying muscle. It's especially beneficial near joints, bony structures, or in areas where soreness or inflammation limits pressure tolerance.

Mid-range strokes (11–13 mm)

Engage both superficial and intermediate tissue layers to increase mobility while preserving comfort. This depth is ideal for most recovery needs—stimulating blood flow, supporting lymphatic drainage, and maintaining pliability in medium-density muscles.

Longer strokes (14–16 mm)

Drive percussion deeper into fascia and thick muscle fibers, reaching embedded tension and fibrotic buildup. This facilitates myofascial release, restores range of motion, and accelerates post-exertional recovery by prompting a deeper parasympathetic response in the tissue.

Does stroke length matter with different modes?

The combination of stroke length and massage mode defines both the intensity and intent of your treatment. Each pairing shapes the way your tissue responds:

Short Stroke + Rolling Wave Mode

Rhythmic, low-impact stimulation. Ideal for sensitive areas or general relaxation, this pairing enhances circulation without overwhelming the tissue, making it perfect for warm-ups or end-of-day wind-downs.

Mid-Range Stroke + Falling Mode

A balanced and restorative approach. Use this combo when treating medium-density muscle groups like hamstrings or delts. It promotes relaxation while gradually tapering off intensity to help the nervous system return to a state of equilibrium.

Long Stroke + Rising Mode

Deep, escalating pressure that mimics progressive sports massage. Perfect for larger or smaller muscles, such as the glutes or quads, this combination facilitates thorough flushing, tension release, and a profound parasympathetic shift post-workout.

In short: depth matters—and learning to match it with the right rhythm unlocks the full potential of Flowgun Ultra.

4. Arm & Positioning

Flowgun Ultra's arm can rotate up to for improved ergonomics. You can lock the arm into different angles to better access hard-to-reach areas.

Try rotating the arm when working on:

- **Lower Back** - So you can reach across the body without straining your shoulder.
- **Hamstrings and Calves** - Adjust for better leverage when seated.
- **Upper Traps and Mid-Back** - Easier access while standing or lying down.

Comfort and angle matter. Adjusting the arm can help reduce tension in your own body while you treat another area.

5. Choosing the Right Head

Each head attachment is engineered for a specific purpose. Using the right one amplifies your results and minimizes unnecessary strain.

Heating Head

Combines percussion with targeted heat up to 50°C. Activates automatically during use to increase circulation and soften stiff muscles. Ideal for use on larger muscle groups. When the ring is red, the heating function is on.

5. Choosing the Right Head

Trigger Cone

Designed for deep, focused release in compact areas like glutes, feet, back, and forearms. Ideal for targeting trigger points and relieving built-up tension.

Fork Head

A dual-point head is designed to stimulate fascia lines with precision and spacing. Ideal for treating larger muscle areas where dual contact points help engage the fascia more effectively.

Axe Head

Angled and flat-edged, this head mimics scraping or flushing techniques across the IT bands, shoulders, back, and lateral muscle groups.

Silicone Hat

A soft, impact-absorbing surface that cushions each percussion. Ideal for use on sore, tender, or sensitive muscle groups — or entry-level treatment.

Flat Head

Broad and rigid, suited for general muscle groups such as thighs, calves, and glutes. Delivers even pressure without overpowering sensitive areas.

Cup Head

Rounded and evenly pressured, perfect for joint-adjacent zones and safe full-body relief.

6. Best Practices & Safety

- Start with low amplitude and speed until your body adapts.
- Do not use over joints, bones, or the cervical spine.
- Avoid prolonged use (more than 3 minutes) on the same spot.
- Clean your massage heads with a damp cloth after each session. Do not submerge the device.

7. Storage & Charging

- Charge the device fully before first use (approx. 3 hours).
- Use the included USB-C charger and adapter. Charging with incorrect voltage may trigger a fault.
- A full charge provides up to 180 minutes of use, depending on the speed and pressure settings.
- Store the Flowgun Ultra and accessories in the included case in a dry environment.

8. Summary

Flowgun Ultra is not just a massage gun—it's a precision recovery tool. By learning its modes, mastering its configurations, and selecting the right tools for your needs, you unlock professional-grade treatment at home or on the move.

Use it purposefully. Adjust it intelligently. And let each session bring you closer to restored, resilient muscles.

For Muscle Soreness

Flowgun Ultra is built to relieve—not just mask—muscle soreness. Whether you're dealing with post-workout stiffness, tension from long hours at a desk, or residual fatigue from high-volume training, the right approach with this device can make all the difference.

Choose the Right Stroke Length

- **Shorter Amplitudes (8–10mm)**
Use when working over sensitive or bony areas. A shorter stroke means less tissue depth, which is ideal for sore calves, forearms, or upper traps.
- **Mid-Range Amplitudes (11–13mm)**
Perfect for general recovery of larger muscle groups like quads, glutes, or hamstrings without overwhelming the tissue.
- **Longer Amplitudes (14–16mm)**
For deep, accumulated soreness in large, dense muscle areas. This setting penetrates further into the fascia and is best used once the tissue is warmed up.

Adjust the stroke using the dial on the side of the device. Real-time feedback from your body matters—lower if it feels sharp, increase if you're not getting enough pressure.

Select the Ideal Mode

- **M2 (Rising Mode)**
A gradual build-up in speed helps slowly bring blood flow into stiff tissue. This is ideal at the start of your session when muscles are cold and sensitive.
- **M4 (Rolling Wave Mode)**
Alternating speeds simulate manual therapy—ramping up and backing off to reduce irritation and flush metabolic waste. This is effective for deep soreness without overstimulating the tissue.
- **M1 (Manual Control)**
Gives you complete control over intensity. Start at a low gear and increase only as the muscle warms up. This is particularly effective when soreness is isolated to specific areas.

Choose the Right Head

- **Silicone Hat**
For tenderness, bruising, or heightened sensitivity. Its soft material reduces impact while still increasing circulation.
- **Heating Head**
Especially effective when tight, stiff muscle fibers cause soreness. The heat component softens tissue and promotes fluid exchange.
- **Flat Head**
For broader surface coverage. Use this over quads, hamstrings, or lats for evenly distributed percussion that avoids pinpoint discomfort.

Application Tips

- **Stay Dynamic**
Move the gun slowly across the muscle, rather than hovering in one spot. One inch per second is a good pace.
- **Limit Duration**
60–90 seconds per muscle group is enough. Spending too long on a sore area may aggravate it.
- **Work Surrounding Areas**
If your quads are sore, spend time working on the glutes and hip flexors as well. This helps relieve systemic tension.

Appropriately used, Flowgun Ultra doesn't just make soreness tolerable—it accelerates recovery and helps restore balance across muscle chains.

For General Recovery

Flowgun Ultra isn't just for treating soreness—it's a daily tool for keeping your muscles supple, your circulation active, and your nervous system balanced. Used consistently, percussion therapy enhances muscle resilience, supports flexibility, and reduces the likelihood of injury or stiffness resulting from everyday activities.

Set the Stroke Length for Maintenance

- **8–10mm**
Ideal for recovery-focused use. This amplitude stays within superficial fascia layers, gently increasing blood flow without irritating the tissue.
- **11–13mm**
Offers moderate depth. Choose this when working over legs, hips, and upper back—areas that carry daily load and benefit from deeper mechanical stimulation.

Avoid long amplitudes unless you are addressing specific deep-tissue needs. The goal here is gentle mobilization, not aggressive treatment.

Use Low-to-Moderate Modes

- **M3 (Falling Mode)**
Begins with higher speeds and gradually slows down. This is great for end-of-day use—starting strong to stimulate, then tapering off to signal the nervous system to relax.
- **M4 (Rolling Wave Mode)**
The alternating rhythm helps stimulate circulation and lymphatic flow while keeping the experience comfortable and balanced.
- **M1 (Manual Mode)**
Choose this when you want precise control over intensity and speed. Gears 2–4 are typically ideal for general recovery.

Best Head Attachments for General Use


- **Cup Head**
Designed for all-around application. The flexible rim conforms to the body and feels intuitive over most areas, including the shoulders, back, legs, and arms.
- **Silicone Hat**
If you prefer a softer touch, this head is ideal for winding down after a long day or for gentle movement or stretching.
- **Flat Head**
For larger muscle groups, such as the thighs and glutes. It delivers an even impact across wide surfaces.

Recovery Timing and Best Practices

- **Use Post-Activity or in the Evening**
This supports parasympathetic activation, preparing the body for rest and regeneration.
- **Session Duration**
Sixty seconds per area is typically sufficient. A full-body sweep might take 8–12 minutes.
- **Stay Consistent**
Building cumulative recovery benefits over time requires three to five sessions per week.

Think of Flowgun Ultra as part of your ongoing care, like brushing your teeth or foam rolling. It keeps tissues hydrated, muscles loose, and your body more adaptable to the physical demands of everyday life.

Before, During, and After Your Workout

Flowgun Ultra is a performance tool, not just a recovery device. Used correctly, it prepares the body for movement, sustains performance, and accelerates post-training repair. This section outlines how to integrate percussion therapy throughout your training timeline to improve warm-ups, enhance intra-session recovery, and achieve more effective cooldowns.

Before Your Workout

**Activate muscles.
Increase range.
Reduce injury risk.**

Using Flowgun Ultra before exercise helps transition muscles from a resting state to a state of active readiness. It raises tissue temperature, improves mobility, and activates neural pathways, particularly in the hips, glutes, and thoracic spine.

- **Stroke Length**
8–14mm for deep tissue priming.
- **Mode**
M2 (Rising Mode) gently ramps up intensity and prepares the nervous system.
- **Time**
30–60 seconds per muscle group (quads, glutes, calves, shoulders).
- **Heads to Use**
Heating head (for glutes/quads), Axe Head (IT band), Trigger Cone (deep hip rotators).

Pro tip: Target both agonist and antagonist muscle groups—e.g., hamstrings and quads—for better balance and control.

During Your Workout

**Reduce tightness.
Maintain performance.
Prevent early fatigue.**

For longer or high-intensity sessions, percussion therapy can be used during breaks or between sets to keep muscles relaxed, flush lactate, and preserve range of motion.

- **Stroke Length**
8–10mm to stay light and fast.
- **Mode**
M1 (Manual Mode) at low gears (2–3) for short bursts of targeted relief.
- **Time**
15–30 seconds on any muscle that starts feeling stiff.
- **Heads to Use**
Cup Head (universal), Silicone Hat (for sensitive or tender zones).

Avoid overuse—this is a tactical reset, not a second workout.

After Your Workout

Flush metabolic waste. Decrease DOMS. Start healing.

Post-exercise is where Flowgun Ultra shines. Percussion therapy stimulates lymphatic flow, eases muscular tension, and signals the parasympathetic system to begin recovery.

- **Stroke Length**
12–16mm to reach deep muscle layers.
- **Mode**
M3 (Falling Mode) or M4 (Rolling Wave Mode) helps shift the body into recovery mode.
- **Time**
60–90 seconds per group, longer if stiffness is present.
- **Heads to Use**
Heating Head (to reduce stiffness and encourage circulation), Fork Head (aductors), Flat Head (for glutes/thighs).


Use slower passes, lighter pressure, and focus on breath to maximize the restorative effect.

hello@flowlife.com
www.flowlife.com
Luntmakargatan 90
113 51 Stockholm

Manual

Flowgun Ultra

1. Rotating Arm
2. Interchangeable Massage Heads
3. Variable Amplitude Wheel
4. Arm Release Button
5. Control Display
6. USB-C Charge Input

ENGLISH

SAFE USE GUIDELINES

- Use only on muscle tissue. Avoid the head, neck, joints, spine, or other bony areas
- Do not apply for more than 60 seconds on one area. Move the device frequently
- Begin with a lower speed and amplitude. Increase gradually based on comfort
- Stop immediately if pain, irritation, or adverse effects occur
- Do not use on injuries, bruises, wounds, or rashes
- Keep hair, loose clothing, and jewelry away during use
- Do not use while sleeping, bathing, or near water sources
- Not suitable for unsupervised children or those with medical restrictions
- Keep vents clear. Never operate on soft or covered surfaces
- Turn off the product before changing heads, cleaning, or charging
- Allow 30 minutes of cooldown after 60 minutes of use
- If the overheating icon appears on the screen, the battery is too hot. The device will automatically activate a safety lock. Let it rest for at least 30 minutes before using again
- Only use approved Flowlife accessories and a charger

CHARGING & POWER

- Charge fully before first use using the included USB-C charger
- Charging time: ~3 hours
- Battery life: Up to 180 minutes (depends on speed and pressure)
- Do not charge overnight or with unapproved adapters (risk of fault)
- Charging indicators are visible on the OLED display

MAINTENANCE & CLEANING

- Wipe with a damp cloth after use. Avoid alcohol-based cleaners
- Do not immerse the device or heads in water
- Store in a cool, dry place using the included case
- If not used regularly, recharge the battery monthly to maintain optimal health

BATTERY & DISPOSAL

- Contains an internal Li-ion battery
- Recycle through certified e-waste channels only
- Do not puncture, short-circuit, or expose to high heat
- Keep out of reach of children

RISK GROUPS – CONSULT YOUR DOCTOR FIRST

- Users with pacemakers, heart disease, or metal implants
- Individuals with epilepsy, nerve disorders, or migraines
- Pregnant individuals
- Recent surgery, injury, thrombosis, or spinal conditions
- If on medication (refer to your prescription leaflet)

INSTRUCTIONS FOR USE

- 1. Start / Stop**
 - Hold the power button (2 sec) to turn on/off
 - The OLED screen will display speed level, mode, and battery status
- 2. Speed Control**
 - Use the Speed (+ / -) buttons to toggle through 6 levels:
1500 → 1700 → 1900 → 2100 → 2300 → 2500 RPM
- 3. Switch Modes**
 - Hold the “M” button to enter mode selection
 - Tap “M” again to cycle through 4 modes:
 - M1: Manual Mode
 - M2: Rising Mode
 - M3: Falling Mode
 - M4: Rolling Wave Mode
- 4. Adjust Stroke Depth (Amplitude)**
 - Use the dial on the rotating arm to set stroke depth from 8–16 mm:
 - 8–10 mm: light, surface-level stimulation (warm-up, sensitive areas)
 - 11–13 mm: balanced recovery (back, hamstrings)
 - 14–16 mm: deep release (glutes, quads, dense tissue)

5. **Arm Positioning**
 - Press and rotate the adjustable arm (up to 90°) for better reach
 - Ideal for accessing the back, hamstrings, and shoulders without strain
6. **Changing Heads**
 - Turn off the device
 - Gently pull to remove the head, and press firmly to insert a new one
 - Heating Head auto-activates when attached and powered on
 - Heat max: 50°C, supports muscle softening and circulation
7. **Included Attachments**
 - Heating Head
 - Trigger Cone
 - Fork Head
 - Axe Head
 - Flat Head
 - Cup Head
 - Silicone Hat

WARRANTY (2 YEARS)

- Covers manufacturing defects
- Excludes: misuse, accidents, unauthorized repairs
- Use outside Flowlife instructions
- Lack of proof of purchase
- Damage exceeding product value

TECHNICAL SPECIFICATIONS

- Amplitude: 8–16 mm (adjustable)
- Vibration Frequency: ≤60 Hz
- Speed Levels:
 - M1 (Manual Mode): 6 adjustable levels
 - RPM range: 1500 – 2500 RPM
- Operation Modes:
 - M1: Manual Mode
 - M2: Rising Mode (Auto Mode)
 - M3: Falling Mode (Auto Mode)
 - M4: Rolling Wave Mode (Auto Mode)
- Stroke Timing: 10-minute intelligent session timer
- Rated Power: 25 W
- Rated Voltage: 21.9 V
- Battery Type: Rechargeable lithium-ion
- Battery Capacity: 2600 mAh
- Battery Runtime: Up to 1–3 hours (varies by mode/force/intensity)
- Charging Time: Approx. 3 hours
- Charging Port: USB Type-C
- Adapter Input: 100–240V ~ 50/60Hz, 0.5A
- Adapter Output:
 - 5V □ 3A
 - 9V □ 2.23A
 - 12V □ 1.67A
- Display: OLED screen showing speed, mode, battery level
- Arm Adjustability: Elbow arm rotation 0–90°
- Auto-Off Safety:
 - 10-minute max operation per session
 - 60 seconds auto-standby when inactive
 - Overheat icon triggers 30-minute cooldown if device gets too hot
- Dimensions:
 - Weight: 1.4kg

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
A close-up, low-angle shot of a person running on a grey, textured surface. The person is wearing a black sports bra and black shorts. Their right leg is extended forward, and their right hand is visible, reaching down. The background is blurred, suggesting motion. The text "Active Lifestyle." is overlaid in white, sans-serif font in the center of the image.

Active Lifestyle.

A blurred photograph of a person running on a paved path. The person is wearing a black sports top and a black wristband. The background is a blurred grey pavement.

Made Possible.

Follow us on Instagram for more inspiration.

A black and white photograph showing a close-up of a person's neck and shoulder. The person is wearing a dark, possibly black, strap over their shoulder. The skin is smooth and has a slight sheen. The word "FLOWLIFE" is overlaid in white, uppercase, sans-serif font across the center of the image.

FLOWLIFE