Flowtens Feet

FLOWLIFE

Brochure

Flowtens Feet

TENS Foot Massager

Table of Content

4 —	Manifest
6 —	Why Flowtens Feet
7 —	Accessories
8 —	Flowtens Feet - A Comprehensive Guide
10 —	Starting your Journey
14 —	For Muscle Soreness
16 —	For General Recovery
18 —	Before, During, and After your Workout
23 —	Contact Information

Manual

- English
- Svenska
- Danska
- Norsk
- Suomi
- Deutsch
- Polski
- Français
- Italiano
- Español
- Nederlands

Manifest

In a world where limits are meant to be shattered, we fuel the fire of relentless pioneers. Our brand champions those who break the mold and embrace the extraordinary, symbolizing a life lived without compromise.

Committed to relentless innovation, we design for those who demand excellence and performance, understanding that the journey into the unknown often leads to unforgettable achievements.

We want to power dreams, aspirations, and triumphs of the audacious, sparking the passions of the unstoppable.


For those who redefine boundaries and explore new horizons, we are Flowlife. Made Possible.

Why Flowtens Feet?

- **Clinically Proven Stimulation, Adapted for Daily Life** - Flowtens Feet uses TENS (Transcutaneous Electrical Nerve Stimulation) to block pain signals and calm irritated nerve endings. It also delivers EMS (Electrical Muscle Stimulation) to engage deep foot muscles and improve postural stability and blood flow. These are trusted technologies in physiotherapy—now adapted for home recovery, wherever you are.
- **Medically Certified, Professionally Designed** - Flowtens Feet is a Class II medical device, tested and certified for safety and effectiveness. It delivers therapeutic-grade stimulation through the soles, helping manage pain, fatigue, swelling, and muscle stiffness—without requiring any movement or manual effort.
- **Recovery Without Resistance** - Using Flowtens Feet requires no pressure, no positioning, and no routine to learn. Moisten your feet, sit down, and let the stimulation do the work. The electrical impulses are delivered passively while you rest, making it ideal for people who need consistent relief but don't have the energy or time to stretch, roll, or mobilise.
- **10 Modes. 50 Levels. One System** - Choose from 10 structured stimulation programs, each with a distinct therapeutic pattern—from acupuncture-style pulses to deep kneading rhythms. Adjust across 50 intensity levels to suit your sensitivity, goals, or time of day. You're always in control—and you can change the feel of your session at any time.
- **Ground-Up Relief** - Foot fatigue is rarely just about the feet. It affects your calves, posture, sleep, and recovery. Flowtens Feet helps you treat the origin, not just the symptoms, by restoring balance where your body meets the ground.


Accessories

When you receive Flowtens Feet, you're not just getting a device—you're getting a complete, clinically grounded recovery tool that's ready to use straight out of the box

Flowtens Feet stimulation device

A compact, all-in-one unit delivering TENS and EMS stimulation through the soles of the feet. Built with an integrated control panel and LCD screen for selecting modes, adjusting intensity, and tracking session time—no remote, no app, no extra setup.

USB-C Charging Cable

For safe and efficient charging. The device runs on a built-in rechargeable battery—no need for disposable batteries or replacement units.

User manual

Includes setup guidance, mode descriptions, intensity recommendations, safety instructions, and usage suggestions for pain relief, circulation, and postural reset.

Flowtens Feet: A Comprehensive Guide

Flowtens Feet is designed to make recovery easy, repeatable, and safe—whether you're addressing discomfort or staying ahead of it. By combining clinically grounded TENS and EMS technology with a user-friendly interface, it delivers consistent stimulation through the soles while you rest. No cables, no effort—just structured relief, one session at a time.

Getting Started: Preparing for Your First Session

Moisten your feet


To ensure proper conductivity, the soles of your feet must be slightly damp before use.

- Use clean water with a spray bottle, a damp towel, or your hands
- Focus on the entire sole, heel, arch, and forefoot
- Do not use lotion, socks, or dry skin
- Without moisture, the stimulation will not transmit effectively

Position the device

- Place the Flowtens Feet on a hard, flat surface (e.g., floor, mat)
- Sit in a stable chair, knees at 90°, feet fully resting on the stimulation zones.
- Keep both feet in complete contact—do not shift, lean, or stand during use

Powering On and Adjusting Settings

- Press the Power button to activate the device
- Use the “+” button to increase intensity gradually
- Start low and raise only to the point of an active but comfortable sensation
- Use the “M” button to cycle through the 10 stimulation modes
- The LCD screen will display your mode number, stimulation level, and session time

Customizing Your Treatment

Flowtens Feet includes 10 built-in stimulation modes, which can be changed using the “M” button. Each program has a unique rhythm and therapeutic purpose. You can switch modes at any time during the session.

Understanding the Ten Modes

Pushing (Mode 1)

Steady, low-frequency pulses that simulate rhythmic acupuncture pressure. Ideal for blocking pain signals in sore, tired feet after prolonged standing or walking.

Soothing (Mode 2)

Wave-like stimulation is designed to gently activate blood flow. Excellent for relieving swelling, heaviness, or cold feet—especially after sedentary periods or long days.

Kneading (Mode 3)

Deep, rolling pulses that mimic therapeutic massage techniques. Helps relax tense tissue and stimulate underlying muscle fibers for improved flexibility and postural balance.

Tapping (Mode 4)

Fast, repetitive pulses resembling light percussive tapping. Effective in relieving stiffness and discomfort by increasing local circulation and nerve responsiveness.

Scrapping (Mode 5)

Repetitive stimulation that mimics scraping motions is often used in traditional therapies. Designed to stimulate the surface nerves and promote relief in tired soles.

Squeezing (Mode 6)

Short, robust EMS bursts that create a compressive squeezing effect. Especially effective for targeting nerve-related discomfort and deep foot tension.

Massage (Mode 7)

Balanced, rhythmic stimulation that helps ease built-up tension in the feet. Ideal for reducing daily stress, calming the nervous system, and restoring a sense of lightness after long periods of standing or walking.

Deep Kneading (Mode 8)	Firm, consistent pressure that mimics manual massage techniques. Designed to ease heavily loaded or strained feet by penetrating deeper muscle groups and connective tissues.
Compression (Mode 9)	Contracting, pulse-driven compression to promote full-foot muscle relaxation. Perfect for resetting after activity or easing lingering tightness in the lower legs and feet.
Acupuncture Squeezing (Mode 10)	A hybrid pattern combining acupuncture-style pulses with deeper activation bursts. Supports neuromuscular reactivation and improves nerve feedback and tone in sensitive feet.


Adjusting Intensity Safely

Flowtens Feet supports 50 intensity levels across all modes. Each level increases the amplitude of the electrical signal, which may feel different depending on the waveform.

- Start at Level 1–5: Suitable for new users, sensitive feet, or nerve calming.
- Progress to 6–30: Ideal for circulation, postural fatigue, or muscle tightness
- Advance to 31–50: For deep tissue recovery, post-training strain, or strong muscle engagement

Important: When switching between modes, intensity does not reset. Always reduce intensity before changing modes, and adjust upward again gradually.

Ending Your Session

- Flowtens Feet automatically shuts off after 30 minutes
- To end earlier, press the Power button manually
- Remove your feet from the pads
- Wipe the surface with a damp cloth (no alcohol or chemicals)
- Allow to air dry before storing
- Charge as needed using the supplied USB-C cable

Consistency builds results

Use Flowtens Feet regularly—1–2 times daily, or a few sessions per week—depending on your needs. The device is most effective when integrated into a comprehensive recovery routine, rather than used solely in response to pain.

For Muscle
Soreness

When soreness builds up in your feet and lower legs, after long shifts, heavy training, or accumulated daily load, Flowtens Feet offers a structured, non-invasive solution. By delivering targeted electrical impulses through the soles, it helps relieve pain, stimulate recovery, and re-engage tired or overloaded foot muscles.

Understanding Foot Muscle Soreness

The feet contain dozens of small stabilising muscles responsible for balance, propulsion, and postural control. These muscles fatigue faster than larger muscle groups and are often overlooked in recovery routines.

Soreness can stem from:

- Prolonged standing or walking
- High-impact training (e.g., running, HIIT, court sports)
- Tight footwear or reduced circulation
- Compensation from postural imbalance is higher up the chain

TENS helps block pain signals, while EMS gently reactivates dormant or fatigued muscle fibres—creating a dual action that relieves discomfort and restores control.

Recommended Modes for Muscle Soreness

Pushing (Mode 1)

Steady, acupuncture-like pulses that interrupt pain signals and calm the nervous system.

- Duration: 20–30 minutes
- Intensity: Start at **5–15**, increase gradually
- Best for: general soreness after standing, walking, or full-body fatigue

Kneading (Mode 3)

Rolling, massage-like stimulation that works deeper into muscular structures.

- Duration: 15–20 minutes
- Intensity: 15–30, depending on tolerance
- Best for: strain from training, overuse, or rigid arches

Compression (Mode 9)

Firm, full-foot contractions that help reset tension and flush out residual fatigue.

- Duration: 15–25 minutes
- Intensity: 20–35
- Best for: end-of-day tightness, heavy muscle load, or reduced circulation

What to Expect

- A noticeable but comfortable tingling or pulsing in the soles
- Gentle twitching or contraction in the calves is normal at higher levels
- Warmth, loosening, and lightness in the feet as the session progresses
- Relief that continues to build after the session ends

Tip: Use Flowtens Feet after long walks, workouts, or travel days to prevent soreness from escalating into dysfunction. When used regularly, it supports healthier loading patterns and faster post-strain recovery.

A black and white photograph of a person's legs bent and crossed, with the text "For General Recovery" overlaid. The person is wearing dark shorts. The background is a plain, light color.

For General
Recovery

Not all discomfort begins with pain. Tightness, heaviness, or postural fatigue often builds up silently. Flowtens Feet helps you address this low-level strain before it escalates by promoting circulation, restoring foot sensitivity, and reactivating dormant muscles. Whether it's the end of a long day or the start of a new one, this is your foundation for full-body recovery.

Why General Recovery Matters

Your feet are under constant demand—even when standing still. Prolonged standing, poor footwear, hard surfaces, or long periods of inactivity can all impair circulation and weaken neuromuscular responsiveness. This leads to:

- Heaviness in the lower legs
- Reduced balance and proprioception
- Sensory fatigue
- Poor recovery between training and workdays

TENS and EMS work together to reverse this cycle—restoring input, enhancing output, and rebalancing activation from the ground up.

Recommended Modes for General Recovery

Soothing (Mode 2)

Wave-like stimulation to gently boost blood flow and reduce pooling in the feet.

- Duration: 10–20 minutes
- Intensity: 5–20
- Best for: morning circulation, post-sitting, or end-of-day reset

Scrapping (Mode 5)

Surface-level nerve activation that enhances foot sensitivity and relieves tiredness.

- Duration: 15–25 minutes
- Intensity: 10–25
- Best for: flat feet, cold feet, or neurological fatigue

Massage (Mode 7)


Rhythmic, balanced impulses that calm the nervous system and ease tension.

- Duration: 15–30 minutes
- Intensity: 10–30
- Best for: everyday decompression, post-work stress, or sensory reset

What to Expect

- A calming, non-intrusive sensation under the soles
- Increased warmth, circulation, and awareness in the feet
- Reduced leg heaviness and postural fatigue
- A subtle feeling of lightness and balance that carries into your next activity

Tip: Use in the morning to improve readiness and stability, or in the evening to unwind tension from the day. Over time, regular sessions help maintain better balance, coordination, and comfort, particularly for individuals who spend extended periods sitting or standing.

Before, During,
and After
Workout

Whether you're lifting, running, cycling, or walking, your feet absorb every rep. Flowtens Feet helps you prepare them before effort, keep them active during recovery phases, and restore them after load. Used consistently, it supports performance, stability, and neuromuscular longevity from the ground up.

Before Workout

Activate foot muscles and reset sensory input

Most warmups focus on large muscle groups, but performance begins with stability. The foot contains small, fast-responding muscles that control balance and ground reaction. Flowtens Feet helps prime these systems before you train.

Benefits

- Stimulates postural and stabilising muscles in the feet
- Increases sensory input and balance awareness
- Improves blood flow, tissue temperature, and oxygenation
- Helps reduce stiffness caused by inactivity (e.g., sitting before training)

Recommended Mode: Kneading (Mode 3)

- Duration: 10–15 minutes seated
- Intensity: 10–25
- Use just before training to increase readiness and lower injury risk

During Workout Recovery

Keep the nervous system engaged between sets

For longer training sessions or circuit-based work, Flowtens Feet can be used passively between blocks to promote recovery without complete rest.

- Maintains light circulation
- Keeps motor units responsive
- Reduces the risk of cold or stiff lower extremities mid-session

Recommended Mode: Scrapping (Mode 5)

- Duration: 5–10 minutes, passive seated
- Intensity: **5–15**
- Best used between sets, supersets, or post-cardio blocks

After Workout

Flush out metabolic waste and restore mobility

After load, tissues are inflamed, fatigued, and often tight. Using Flowtens Feet immediately post-exercise helps reset the neuromuscular system and accelerate recovery.

Benefits

- Stimulates lymphatic flow and venous return
- Reduces residual stiffness and inflammation
- Supports calm-down and parasympathetic recovery
- Helps re-establish proper foot strike and mobility after an effort

Recommended Modes

- Pushing (Mode 1) – for pain relief and calming
- Compression (Mode 9) – for deep decompression
- Duration: 20–30 minutes seated
- Intensity: **15–35**, depending on effort level and tolerance


Tip: Using Flowtens Feet around your training supports long-term performance by reducing tension buildup and improving recovery efficiency. Over time, it helps protect against overuse patterns that begin in the feet and cascade upward into the knees, hips, and spine.

hello@flowlife.com
www.flowlife.com
Luntmakargatan 90
113 51 Stockholm

Manual

Flowtens Feet

1. Foldable Stand
2. USB-C Port
3. LCD Display
4. Feet soles
5. Control Panel
6. Mode Button
7. Minus
8. On/Off Button
9. Plus

ENGLISH

GETTING STARTED SAFELY

To ensure optimal performance and personal safety, please follow the guidelines below when using Flowtens Feet:

- Only use on clean, bare, slightly damp feet
- Do not apply lotions, oils, or wear socks during use
- Use the device while seated
- Never stand or apply full body weight to the device
- Avoid use near water or while bathing
- Flowtens Feet is not water-resistant and should not be used in wet environments
- Keep the device dry
- If it becomes wet, allow it to dry completely before reuse
- Do not use while sleeping, driving, or operating machinery

HEALTH AND MEDICAL PRECAUTIONS

Consult a healthcare provider before use if you:

- Have a pacemaker, defibrillator, or any implanted medical device
- Are pregnant
- Have epilepsy or other neurological disorders
- Suffer from heart disease, arrhythmias, or circulatory disorders
- Are recovering from recent surgery
- Have diabetes with complications (e.g., peripheral neuropathy)
- Are under treatment for chronic conditions such as liver disease or gastrointestinal bleeding
- Experience sensory impairments (e.g., reduced sensitivity to heat, pain, or pressure)
- Have skin conditions such as eczema, dermatitis, or open wounds on the feet
- Not recommended for children without prior medical consultation

SAFE USE GUIDELINES

Press the Power button to activate the device

- Use the “+” button to increase intensity gradually
- Use the “M” button to cycle through the 10 stimulation modes
- The LCD screen will display your mode number, stimulation level, and session time

Session Duration:

- Limit use to 30 minutes per session
- The device will shut off automatically

Frequency:

- Use up to twice daily, with sufficient rest between sessions

Pressure:

- Let the electrical impulses do the work – do not press with excess weight

Intensity Levels:

- Always begin at the lowest setting and increase gradually to avoid discomfort or overstimulation

DISPOSAL

- Flowtens Feet has a battery inside
- Do not dispose of batteries in household waste or water
- Do not puncture, burn, or short-circuit batteries
- Follow local regulations for proper disposal or recycling

WARRANTY & SUPPORT

Flowtens Feet is covered by a 2-year warranty against material and manufacturing defects from the date of purchase

If your device fails within this period, we will repair or replace it at no additional charge

This warranty does not cover:

- Damage resulting from misuse, accidents, or improper handling
- Unauthorized repairs or modifications
- Use that does not comply with this user guide
- Damage from improper storage or environmental exposure
- Lack of proof of purchase

ELECTRICAL SPECIFICATION

- Product size: 330×337.5×42.8 mm
- Rated: 3.7V DC 100mA
- Pulse repetition frequencies: 0–200 Hz
- Maximum amplitude of output voltage: 72V
- Battery capacity: 400mAh
- Rated frequency: 50–60 Hz
- Rated output: 5V 1000mA
- Pulse duration: 120 μ s
- Effect of load impedance: 1k Ω

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WERKOMSTANDIGHEDEN

- 10 – 40°C
- 30 % RH ~ 75 % RH
- 860 hPa tot 1060 hPa

Symbol		Description	
	Manufacturer		Consult instructions for use
	Authorized representative in the European Community		Authorized representative in the UK
	Date of manufacture		Class II symbol
	Batch code		Type BF according to electrical safety requirements
	Serial number		Recycling symbol
	Keep dry		Caution
	No sitting		No stepping on surface
	Operating temperature limitation 0°C to 40°C (32°F to 104°F)		
 2862	Complies with applicable EU Standards		
IP21	Protected from touch by fingers and objects greater than 12 millimeters, and protected from condensation.		
	Indicates the item is a medical device		

Specifications

Model number: AST-300Q

This equipment belongs to Type BF of Class II according to Clause 5.

Shenzhen OSTO Medical Technology Company Limited
2104, Unit4, Building2, Zhongxin Longsheng Square,
NO. 219 Shenshan Road (Longgang Section), Longxin Community,
Baolong Street, Longgang District, 518116 Shenzhen City, Guangdong Province, China

Phoenix Medtech GmbH
Koenigsberger Strasse 11, 64839,
Muenster Hessen, Germany
Contact: Kevin Han
Tel: +49 6071 4977513
Email: contact@phoenix-medtech.de

UK Responsible Person

SHARE INFO LTD
3rd Floor, Office C, Townend House, Park Street,
Walsall, West Midlands, WS1 1NS, UK.

Flowlife.com
Luntmakargatan 90
113 51 Stockholm
hello@flowlife.com